
# PROVOX LaryTube

**Rx  
ONLY**


Figure 1


Figure 2


Figure 3


Figure 4


Figure 5


Figure 6a


Figure 6b


Figure 6c


Figure 6d


Figure 6e


Figure 6f


Figure 6g


Figure 6h

## **Prescription information**

**CAUTION:** United States Federal law restricts this device to sale, distribution and use by or on order of a physician or a licensed practitioner. The availability of this product without prescription outside the United States may vary from country to country.

## **Disclaimer**

Atos Medical offers no warranty—neither expressed nor implied—to the purchaser hereunder as to the lifetime of the product delivered, which may vary with individual use and biological conditions. Furthermore, Atos Medical offers no warranty of merchantability or fitness of the product for any particular purpose.

## Contents

EN - ENGLISH .....	7
DE - DEUTSCH .....	14
NL - NEDERLANDS .....	22
FR - FRANÇAIS .....	29
IT - ITALIANO .....	36
ES - ESPAÑOL .....	43
PT - PORTUGUÊS .....	50
SV - SVENSKA .....	57
DA - DANSK .....	64
NO - NORSK .....	71
FI - SUOMI .....	78
IS - ÍSLENSKA .....	85
ET - EESTI .....	92
LV - LATVIEŠU .....	99
LT - LIETUVIŲ KALBA .....	106
CS - ČESKY .....	113
HU - MAGYAR .....	120
SL - SLOVENŠČINA .....	127
SK - SLOVENČINA .....	134
PL - POLSKI .....	141
RO - ROMÂNĂ .....	148
HR - HRVATSKI .....	155
SR - SRPSKI .....	162
EL - ΕΛΛΗΝΙΚΑ .....	169
TR - TÜRKÇE .....	177
KA - ქართული .....	184
HY - ՀԱՅԵՐԷՆ .....	193
RU - РУССКИЙ .....	201
ID - BAHASA INDONESIA .....	209
MS - BAHASA MELAYU .....	217
JA - 日本語 .....	224
KO - 한국어 .....	231
ZHTW - 繁體中文 .....	238
TH - ไทย .....	244
VI - TIẾNG VIỆT .....	251
HE - עברית .....	263
AR - العربية .....	270
FA - فارسی .....	277

## 1. Intended Use

Provox LaryTube is a holder for devices in the Provox HME System intended for vocal and pulmonary rehabilitation after total laryngectomy.

For patients with a shrinking tracheostoma, it is also used to maintain the tracheostoma for breathing.

Provox LaryTube is intended for single patient use.

## 2. Contraindications

Provox LaryTube is not intended to be used by patients who

- are under any form of mechanical ventilation.
- have damaged tracheal or tracheostoma tissue.

## 3. Description of the Device

Provox LaryTube is a tube made of medical-grade silicone rubber. The purpose of the device is to create a comfortable and airtight fit between Provox LaryTube and the tracheostoma, and also to provide attachment for devices from the Provox HME System.

They are delivered single packed, non-sterile.

There are three models: Standard, Fenestrated, and with Ring (see Figure 1).

The standard version and ring version can be fenestrated so that air can go through the voice prosthesis for voice prosthesis users. The holes are punched by using a biopsy punch. The fenestrations should be in a position that will optimize the airflow through the voice prosthesis.

<p><b>Standard versions</b> Made for use with or without a voice prosthesis. Can be attached with a Provox TubeHolder or Provox LaryClip.</p>	See Figure 1a
<p><b>Fenestrated versions</b> For voice prosthesis users. Can be attached with a Provox TubeHolder or Provox LaryClip.</p>	See Figure 1b
<p><b>Ring versions</b> Made for use with or without a voice prosthesis. Can only be attached with a Provox Adhesive.</p>	See Figure 1c

## Provox LaryTube Parts

Tube	See Figure 1d
Shield (conical)	See Figure 1e
HME and accessory holder	See Figure 1f
Wings	See Figure 1g
Ring for adhesive connection	See Figure 1h
Fenestration holes	See Figure 1i

## 4. Warnings

- Do not reuse a used Provox LaryTube in another patient. The device is intended for single patient use only. Reuse in another patient may cause cross-contamination.
- Do not use lubricating gels if the patient has allergies related to these substances.
- Do ensure that the patient has been trained in the use of the device. The patient should demonstrate the ability to understand and consistently follow the Instructions for Use without clinician supervision.

## 5. Precautions

- Always use a Provox LaryTube of an appropriate size. If Provox LaryTube is too wide, too narrow, too long, or too short, it may cause tissue damage, bleeding, or irritation. In addition, it may be harder to produce voice, because Provox LaryTube may obstruct the voice prosthesis. Inserting, removing, and pressing on the HME Cassette may exert pressure on the voice prosthesis if Provox LaryTube is not sized appropriately. A Provox LaryTube that is too narrow may cause the stoma to shrink.
- Do not use force during insertion. Always insert and remove Provox LaryTube according to the instructions below. If the product is pushed too far into trachea, it may obstruct breathing. The use of force during insertion may cause tissue damage, bleeding or irritation, and accidental dislodgement of the voice prosthesis and/or Provox XtraFlange if used.
- If fenestration is necessary, make sure that it is done with Provox Fenestration Punch. This will ensure that the holes are small and smooth. Fenestration holes that are too large may cause dislodgement of the voice prosthesis or Provox XtraFlange. Fenestration holes with sharp edges may cause formation of granulation tissue. Make sure that no silicone residue is left in the holes or lumen of Provox LaryTube after fenestration.

- Do not use a dirty or contaminated Provox LaryTube. Clean and disinfect the device according to the cleaning and disinfection instructions below.
- Do not use a damaged device since this might cause aspiration of small parts or cause formation of granulation tissue.
- Use only water-soluble lubricants. Oil-based lubricants (e.g. Vaseline) may weaken, damage, or destroy the product and should not be used.
- Carefully monitor the tracheostoma tissue during radiation therapy. Stop using the device if the stoma gets irritated or starts to bleed.
- Patients with bleeding disorders or patients on anticoagulants should not use the device if it causes bleeding that reoccurs.
- Stop using the device if granulation tissue forms in the tracheostoma.
- Do not clean the device while inside the stoma since this could cause tissue damage. Always remove the device from the stoma before cleaning.

## 6. How to Use

### 6.1 Preparation

#### Choosing the Correct Size

##### **Diameter:**

Provox LaryTubes are available in four different diameters: 8, 9, 10, and 12. Measure the size of the tracheostoma (largest diameter) with a ruler and fit Provox LaryTube accordingly. In irregularly shaped tracheostomas, it may be difficult to achieve an airtight fit.

##### **Length:**

Provox LaryTubes are available in three different lengths: 27, 36, and 55 mm. To help choose the size(s) of Provox LaryTube to prescribe to the patient, the clinician can use a Provox LaryTube Sizer Kit. The sizer kit contains samples ("sizers") of commercially available Provox LaryTubes.

**Caution:** During fitting, always make sure that the device does not catch on to the tracheal flange of the voice prosthesis (see Figure 3).

#### Preparing for Insertion

Before first use, check that you have the correct size and that the package is not damaged or opened. Do not use otherwise.

- Make sure that your hands are thoroughly clean and dry before placing them in the area of the stoma or handling Provox LaryTube and/or accessories.
- Have a mirror and sufficient light to illuminate the stoma.
- Always inspect Provox LaryTube thoroughly before each use (i.e. before insertion) and ensure that it is undamaged and shows no signs of tear, cracks, or crusts. If you detect any damage, do not use the product and obtain a replacement.

- Before insertion, rinse the device with drinking water to remove any disinfectant residue (see Figure 2). Make sure that you only insert clean devices. In a hospital setting, it is better to use sterile water rather than drinking water.
- If needed, lightly lubricate the tube with a small amount of water-soluble lubricant.

#### **Recommended water-soluble lubricants**

KY Jelly®
Surgilube®
Xylocaine® (lidocaine)

**Caution:** Use only water-soluble lubricants. Oil-based lubricants (e.g. Vaseline) may weaken, damage, or destroy the product and should not be used at any time. Make sure that you do not lubricate the HME holder, HME Cassette, or any accessory which is to be held by Provox LaryTube. This may lead to accidental detachment.

## **6.2 Operating Instruction**

### **Insertion**

#### **Provox LaryTube Standard and Fenestrated**

1. Gently insert Provox LaryTube directly into the tracheostoma (see Figure 3).
2. Attach these tubes with either Provox LaryClip or Provox TubeHolder (see Figures 4b and 4c).

#### **Provox LaryTube with Ring**

1. Attach a Provox Adhesive and then gently insert Provox LaryTube through the adhesive into the tracheostoma (see Figure 4a).
2. Make sure that the adhesive is centered properly so that Provox LaryTube does not press on the sides of the stoma.
3. Insert the Provox HME System component into the HME holder.

**Note:** Some patients might experience more coughing after inserting Provox LaryTube for the first time or when occluding the HME to speak. This usually diminishes as the user adjusts to the device. If discomfort persists, Provox LaryTube should be removed.

## Removing the Device

1. Hold Provox LaryTube in place and remove the Provox HME System component from the tube (see Figure 5).
2. Remove Provox LaryTube from the stoma.

**Note:** Insertion and/or removal of the device may occasionally cause some slight bleeding, irritation, or coughing. Inform your clinician if these symptoms persist.

## 6.3 Cleaning and Disinfection

**Caution:** Do not clean the device while inside the stoma. This can cause injury. Always remove the device from the stoma before cleaning.

Clean Provox LaryTube after each use. Do not use water other than drinking water to clean and rinse the device.

**Caution:** Be careful not to burn your fingers in the warm cleaning water.

1. Rinse Provox LaryTube in drinking water (see Figure 6a).
2. Place the device in non-oil based dish soap (see Figure 6b).
3. Clean the inside of the tube with Provox TubeBrush (see Figures 6c and 6h).
4. If a fenestrated tube is used, the holes should be cleaned with a Provox Brush (see Figures 6d and 6h).
5. Rinse Provox LaryTube in drinking water (see Figure 6e).
6. Disinfect the device at least once a day (see Figure 6f) with one of the following methods:
  - Ethanol 70% for 10 min
  - Isopropyl alcohol 70% for 10 min
  - Hydrogen peroxide 3% for 60 min

**Caution:** Do not use the device until it is completely dry (see Figure 6g). Inhalation of disinfectant fumes can cause severe coughing and irritation of the airways.

If Provox LaryTube looks dirty or has air dried in an area with a risk of contamination, the device should be both cleaned and disinfected before use.

A risk of contamination could exist if the device has been dropped on the floor or if it has been in contact with a pet, someone with a respiratory infection, or any other gross contamination.

During hospitalisation, it is important to both clean and disinfect Provox LaryTube after use but also before use since there is an increased risk of infection and contamination. In a hospital, it is better to use sterile water for cleaning and rinsing rather than drinking water.

**Caution:** Do not clean or disinfect by any other method than described above as this might cause product damage and patient injury.

## 7. Compatible Products

**Caution:** Use only original components. Other accessories may cause product damage and/or malfunction.

### Products compatible with Provox LaryTube

<b>Provox HME Cassettes</b>	Heat and Moisture Exchangers
<b>Provox FreeHands HME</b>	For speaking without manual stoma occlusion.
<b>Provox ShowerAid</b>	Can help to prevent water from entering Provox LaryTube during showering.
<b>Provox TubeHolder or Provox LaryClips</b>	Can be used for additional support to prevent Provox LaryTube from falling out of the stoma.
<b>Provox Adhesive</b>	Is used as attachment and keeping Provox LaryTube with Ring in place.

## 8. Troubleshooting Information

### Bleeding of the Stoma

Stop using the device. After the tracheostoma is healed, gradually try to increase the duration of time when using the device. If even short periods of time cause bleeding or if the bleeding persists after removal of the device, stop using it and consult your clinician.

### Rapid Shrinking of the Stoma

If you are wearing Provox LaryTube because of a rapidly shrinking tracheostoma, discuss with your clinician the possibility of having two devices (of the same size). This would allow you to remove Provox LaryTube from your stoma and replace it immediately with a clean one. Cleaning of the removed device can then be done after you have inserted the clean Provox LaryTube.

### Stoma Enlarges During Use

If your tracheostoma enlarges after using Provox LaryTube, ask your clinician if it is possible for you to use two devices (one with a smaller diameter and one with a larger diameter). This allows you to use the larger size Provox LaryTube when the tracheostoma has enlarged, and the smaller size to allow the tracheostoma to shrink. It may be necessary to discontinue use of the tube during the night so that the tracheostoma can shrink. In that case, discuss alternative options to attach the HME Cassette with your clinician so that you can continue HME Cassette use during the night.

## **9. Storage Instructions**

When not in use, clean and disinfect the device as described above and then store Provox LaryTube in a clean and dry container at room temperature. Protect from direct sunlight.

## **10. Device Lifetime**

Provox LaryTube may be used for a maximum period of 6 months. If the device shows any signs of damage, it shall be replaced earlier.

## **11. Disposal**

Always follow medical practice and national requirements regarding biohazard when disposing of a used medical device.

## **12. Additional Information**

### **12.1 Compatibility with MRI Examination**

#### **MR Safe**

This device does not contain any metallic elements and has no potential for interaction with the MRI field.

## **13. Reporting**

Please note that any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the national authority of the country in which the user and/or patient resides

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Atos and the Atos logo are trademarks of Coloplast A/S.
© 2025 Coloplast A/S. All rights reserved.
Atos Medical AB