# ProTrach<sup>TM</sup>

## DualCare<sup>TM</sup>

### Instructions for Use

# Atos

Atos Medical Your voice


Figure 1

Figure 2a


Figure 2b


Figure 3


Figure 4a


Figure 4b


Figure 5a


Figure 5b


Figure 5c


Figure 6a


Figure 6b


Figure 6c



Figure 6d


Figure 7a


Figure 7b


Figure 8


Figure 9


Figure 10


Figure 11

Figure 12a

Figure 12b


Figure 12c


Figure 12d

Figure 13a

Figure 13b

Figure 13c


Figure 13d


Figure 13e

Figure 14a


Figure 14b

Figure 14c


Figure 14d

Figure 15a


Figure 15b


Figure 16

## **Prescription information**

**CAUTION:** United States Federal law restricts this device to sale, distribution and use by or on order of a physician or a licensed practitioner. The availability of this product without prescription outside the United States may vary from country to country.

## **Disclaimer**

Atos Medical offers no warranty – neither expressed nor implied – to the purchaser hereunder as to the lifetime of the product delivered, which may vary with individual use and biological conditions. Furthermore, Atos Medical offers no warranty of merchantability or fitness of the product for any particular purpose.

## **Patent and trademarks**

ProTrach<sup>®</sup> is a registered trademark and DualCare<sup>™</sup> is a trademark of Atos Medical AB. The products are protected by US patent no. US 6,921,417, JP patent no. 4,156,832 as well as other patents.

# Contents

ENGLISH .....	9
DEUTSCH .....	24
NEDERLANDS .....	41
FRANÇAIS .....	58
ITALIANO .....	75
ESPAÑOL .....	92
PORTUGUÊS .....	109
SVENSKA .....	127
DANSK .....	142
NORSK .....	158
SUOMI .....	173
ORDERING INFORMATION .....	189

The Instructions for Use, which accompanies this product, may be revised from time to time and must therefore be reviewed prior to each procedure in which the product is used.

## Contents

1. Descriptive information.....	10
1.1 Intended use.....	10
1.2 CONTRAINDICATIONS .....	10
1.3 Description of the devices .....	11
1.4 WARNINGS.....	14
1.5 PRECAUTIONS .....	16
2. Instructions for use.....	17
2.1 Preparation and assembly of the speaking devices.....	17
2.2 Insertion and usage of Speaking Valve with HME 15 / 22 .....	18
2.3 Insertion and usage of HME DigiTop and HME DigiTop O <sub>2</sub> .....	19
2.4 Removal of the HME from the speaking devices.....	20
2.5 Pediatric Use .....	20
2.6 Cleaning and disinfection .....	21
2.6.1 Cleaning of the speaking devices.....	21
2.6.2 Disinfection of the speaking devices.....	22
2.7 Storage instructions.....	23
2.8 Device lifetime .....	23
2.9 Disposal .....	23
2.10 Accessories .....	23
3. Additional information .....	23
3.1 Compatibility with MRI Examination.....	23
3.2 Ordering information .....	23
3.3 User assistance information.....	23

## Symbol explanation

Keep away from direct sunlight and heat.

Keep dry and away from rain and other moisture.

# 1. Descriptive information

## 1.1 Intended use

ProTrach DualCare is a combined Speaking Valve and Heat and Moisture Exchanger (HME) intended for spontaneously breathing tracheostomized patients using a tracheostomy tube with a deflated cuff, or a tracheostomy tube without cuff. In HME-mode the device conditions inhaled air by retaining heat and moisture from the exhaled air.

By turning the lid of the Speaking Valve into speaking mode air is re-directed to enable speech.

The entire device is for single patient use and the HME-part is for single use.

**Patient Population:** For spontaneously breathing tracheostomized patients (adults and pediatric patients greater than 10 kg in weight) using a tracheostomy tube with a deflated cuff, or a tracheostomy tube without cuff.

**Environment of Use:** Hospitals, ICU, sub-acute care institutions and home.

## 1.2 CONTRAINDICATIONS

### General

ProTrach HME 15 or 22 in combination with either ProTrach DualCare Speaking Valve or ProTrach HME DigiTop are contraindicated for:

- Use in combination with an in-line ventilator.
- Patients who are unable to handle or remove the device themselves when needed, unless the patient is under constant supervision of a clinician or a trained caregiver. For example, patients who are unable to move their arms, patients with decreased levels of consciousness, or patients with diseases that put them at risk for unpredictable periodic loss of consciousness.
- Patients who cannot tolerate the added dead space of 5 ml, or who cannot tolerate the added breathing resistance of 170 Pa / 1.7 cm H<sub>2</sub>O. This should be evaluated by a clinician.
- Unresponsive or sedated patients. The patient must be responsive and attempting to communicate in order to use the device. The patient should be able to follow instructions.

The HME 15 or 22 in combination with either Speaking Valve or HME DigiTop must NOT be used on a single lumen tube (tube without an inner tube), unless the patient or caregiver is able to reinsert the tube themselves after emergency removal.

### **Speaking Valve specific**

The use of **Speaking Valve** (in combination with HME 15 or 22) is **additionally** contraindicated for the following patient groups:

- Laryngectomized patients since the device will prevent the ability to exhale if the Speaking Valve is unintentionally set to speaking mode.
- Patients suffering from severe aspiration.
- Patients with severe obstruction in the area of the tracheostomy tube or in the upper airways, such as tracheal and/or laryngeal stenosis, since this may cause air trapping.
- Patients with very thick and copious secretions which might block the device.

**DO NOT** use the Speaking Valve:

- In combination with a tracheostomy tube with the cuff inflated. The cuff must be completely deflated before placing and during all use of the Speaking Valve.
- In combination with a tracheostomy tube with a foam cuff.
- In combination with a tracheostomy tube with a self-inflating cuff.
- When the size of the tracheostomy tube does not allow for airflow through the upper airways.
- In combination with an endotracheal tube.

**Use of the Speaking Valve in these circumstances can restrict exhalation through the upper airways and cause suffocation!**

**DO NOT** use the Speaking Valve during sleep since the airway could be blocked unintentionally. During sleep the HME DigiTop (in combination with HME 15 or 22) should be used instead (Figure 2a or 2b).

## **1.3 Description of the devices**

The Speaking Valve, the HME DigiTop, and the HME DigiTop O<sub>2</sub> are all used to enable speaking. In this Instructions for Use they will therefore be referenced to as “speaking devices”.

### **Speaking Valve and HME 15 / 22**

The reusable Speaking Valve is used with a single use 15 mm or 22 mm Heat and Moisture Exchanger (HME). The HME is placed so that it prevents direct contact between the Speaking

Valve and the airways. This prevents the Speaking Valve from being clogged or soiled by mucus.

The Speaking Valve has two modes: speaking mode and HME mode.

In speaking mode, a flexible membrane is positioned in the airflow openings and acts as a one way valve. It opens during inhalation so the patient can inhale through the device. During exhalation, the membrane remains closed and the air is re-directed through the upper airways and the vocal folds. Thereby the patient is able to speak. In speaking mode the inhaled air does not get conditioned since the exhaled air goes out through the upper airways.

In HME mode the membrane is moved out of the way of the airflow so that the patient both inhales and exhales through the device. The inhaled air is conditioned by the heat and moisture that is retained from the exhaled air in the impregnated HME media. In HME mode, speaking is not possible.

The device is switched between the modes by rotating the lid of the Speaking Valve until it clicks into the desired position.

The ProTrach HME 15 is available in two versions: the XtraMoist with better humidification properties and the Regular with lower breathing resistance. The XtraMoist should be used when the slightly higher breathing resistance is tolerated. The Regular should be used when lower breathing resistance is required, for example when first getting used to HME use or during physical activity.

## **HME DigiTop and HME DigiTop O<sub>2</sub>**

The HME DigiTop and the HME DigiTop O<sub>2</sub> enable use of the HME without the Speaking Valve, and can manually be occluded to enable speaking.

The HME DigiTop O<sub>2</sub> provides possibility to connect oxygen tubing of 1/8 inch (3.2 mm) diameter to the Oxygen port connector for patients requiring additional Oxygen (Figure 16).

## ProTrach DualCare Set 15 / Set 22 include:

See Figure 3

- (a) 1 pc ProTrach DualCare Speaking Valve
- (b1 or b2) 30 pcs (b1) ProTrach HME 15 Regular or  
(b2) ProTrach HME 22 Regular
- (c) 1 pc ProTrach Connection Strap clip
- (d) 1 pc ProTrach Connection Strap string
- (e) 1 pc ProTrach HME DigiTop
- (f) 1 pc Removal Aid
- 1 pc ProTrach DualCare Instructions for Use

## ProTrach DualCare Speaking Valve / Speaking Valve Blue package include:

See Figure 3

- (a) 1 pc ProTrach DualCare Speaking Valve or Speaking Valve Blue
- (c) 1 pc ProTrach Connection Strap clip
- (d) 1 pc ProTrach Connection Strap string
- (e) 1 pc ProTrach HME DigiTop or DigiTop Blue
- 1 pc ProTrach DualCare Instructions for Use

## Technical data ProTrach DualCare

Speaking Valve with HME 15		Speaking Valve with HME 22	
Pressure drop at 30 l/min in speaking mode, HME 15 Regular:	125 Pa	Pressure drop at 30 l/min in speaking mode:	136 Pa
Pressure drop at 30 l/min in HME mode, HME 15 Regular:	35 Pa	Pressure drop at 30 l/min in HME mode:	48 Pa
Pressure drop at 60 l/min in HME mode, HME 15 Regular:	123 Pa	Pressure drop at 60 l/min in HME mode:	166 Pa
Pressure drop at 30 l/min in speaking mode, HME 15 XtraMoist:	166 Pa		
Pressure drop at 30 l/min in HME mode, HME 15 XtraMoist:	78 Pa		
Pressure drop at 60 l/min in HME mode, HME 15 XtraMoist:	*		
Dead space: (Speaking Valve incl HME)	4.05 ml	Dead space: (Speaking Valve incl HME)	4.62 ml

\* During increased physical activity that requires lower breathing resistance the Regular HME is recommended. Therefore, the XtraMoist is not tested in these conditions.

Oxygen concentration measured at: Tidal volume 1 L, breathing frequency 10/min, O <sub>2</sub> flow 10 L/min	
DigiTop O <sub>2</sub> with HME 15	70%
DigiTop O <sub>2</sub> with HME 22	71%

## 1.4 WARNINGS

### Speaking Valve specific

(Speaking Valve in combination with HME 15 or 22).

#### **WARNING**

If a tracheostomy tube with a cuff is used, **the cuff must be completely deflated** before placing, and during all use of the Speaking Valve. The patient cannot breathe if the cuff is not completely deflated. There should be sufficient airflow around the tracheostomy tube and the deflated cuff. **DO NOT** use the Speaking Valve for patients who use a tracheostomy tube with a self-inflating or foam-filled cuff.

**Use of Warning Labels:** Attach warning labels provided with Speaking Valve to the tracheostomy tube pilot balloon inflation line (Figure 1).

**Monitoring:** Patients with a tracheostomy tube with a cuff using the Speaking Valve must be monitored according to physician direction. Continuous observation is needed to make sure that the patient can exhale through the upper airway.

**Education:** To ensure correct use of the Speaking Valve, the patient, family, caregivers, and all personnel working with the patient must be instructed on how to use Speaking Valve, including contraindications, warnings, warning labels and precautions.

#### **Do**

- Do make sure the patient has the respiratory ability to exhale around the tracheostomy tube and out of the nasal and oral cavities.
- Do make sure the patient has been cleared by the physician for cuff deflation. Cuff deflation is mandatory with the Speaking Valve to allow exhaled air to pass around the tracheostomy tube and upper airways.
- Do make sure the tracheostomy tube is sized to allow for sufficient airflow around the tracheostomy tube to facilitate speech. The cuff on the tube can also generate an obstruction even when deflated and this should be taken into consideration during airway patency assessment.
- Do evaluate, per physician direction, if changing to a smaller tracheostomy tube is indicated to provide sufficient airflow to allow use of the Speaking Valve.

- Do evaluate patients with thick unmanageable secretions that may cause airway obstructions for use of the Speaking Valve.
- Do evaluate patient experiencing difficulty utilizing the Speaking Valve, as the patient may have airway obstructions due to stenosis, tissue mass, tracheomalacia, granulation, vocal cord paralysis in the midline position, secretions, or a tracheostomy tube that is oversized for the patient's trachea. If obstruction is corrected, the patient should be re-evaluated for the Speaking Valve.

## **General**

(HME 15 and 22 in combination with either Speaking Valve or HME DigiTop).

### **Do**

- Do make sure that patients with severe pulmonary diseases, e.g. lung emphysema, severe asthma etc. are examined and evaluated by a clinician to ensure suitability before they start using the device.
- Do make sure that the patient or caregiver always has the ability to remove the device quickly from the tracheostomy tube if needed. Coughing up a lot of mucus may suddenly block or nearly block the tube or HME.
- Do make sure that the patient, caregiver and others understand the closing function and the use of the holes. Do explain that blocking the holes in the Speaking Valve or the HME DigiTop by clothing, for example, may cause difficulty in breathing or suffocation.

### **Do not**

- Do not bathe or swim with the device (Figure 4a or 4b). This could cause water to enter the airways. The speaking devices do not prevent this in any mode.
- Do not use the device if it has been damaged or contaminated.
- Do not remove the foam from the HME. Using the device without the foam increases the risk of aspiration of small device parts, it may also negatively impact valve functionality due to direct mucus contact, and also without the foam the important HME function is lost.
- Do not reuse the HME, e.g. by washing it. By washing the HME the important HME function is lost. Reuse can also increase the risk of infection.

- Do not use any part of the device for more than one patient since this may cause cross contamination. The device is intended for single patient use only.

## **HME 22 specific**

(HME 22 in combination with either Speaking Valve or HME DigiTop).

### **Do**

- If the Speaking Valve comes loose from the HME 22 during the removal of the device, always remove the inner tracheostomy tube or the single lumen tube together with the HME instead of trying to remove the HME from the tube connector.

Attempts to remove the HME 22 from the connector while the tube is still in place is difficult, and the airways may become blocked by mucus.

## **1.5 PRECAUTIONS**

### **Do**

- Do always check before each use that the Speaking Valve rotates properly between the two modes. The Speaking Valve may get stuck in one of the modes due to the presence of mucus. If this occurs, clean the device thoroughly and if the problem persists stop using the device and replace it.
- Do check the device if you experience increased breathing resistance as this might be a sign of blockage by (e.g., by mucus).
- Do consult your clinician if there is a problem in connection with the device.

### **Do not**

- Do not use humidifiers or heated humidified oxygen via a mask over the tracheostoma while using the device since it will cause the HME to become too wet. If oxygen therapy is required, use only non-heated humidified oxygen.
- Do not administer medicated nebulizer treatment over the tracheostoma while using the device. The medication can become deposited in the HME and medications may also adversely affect the Speaking Valve membrane. If the device is inadvertently used during a medicated nebulizer treatment, the device must be removed immediately. The HME must be discarded, and the Speaking Valve and the

HME DigiTop must be cleaned thoroughly to remove all medication residues.

## 2. Instructions for use

The speaking devices (Speaking Valve, HME DigiTop and HME DigiTop O<sub>2</sub>) are reusable. The HME cassettes are disposable and must be exchanged at least every 24 hours.

### 2.1 Preparation and assembly of the speaking devices

**CAUTION:** Use only original accessories with the device. Use of other accessories may cause product malfunction and harm the patient.

When the Speaking Valve in combination with a HME is not advisable, evaluate the possibility of using the ProTrach HME DigiTop in combination with a HME.

Before each use, please check that the speaking device (Speaking Valve, HME DigiTop, or HME DigiTop O<sub>2</sub>) and HME are undamaged and that the Speaking Valve functions as intended, e.g., it clicks into the HME mode and speaking mode. If the device does not function as intended or seems damaged, do not use the device and replace it.

#### Attach the ProTrach Connection Strap

The ProTrach Connection Strap can be used to secure the Speaking Valve to the neckband (optional). This is done to avoid that the device is lost, for example, when coughing or sneezing.

1. Insert the string of the Connection Strap into the clip (Figure 5a and 5b).
2. Insert the string of the Connection Strap to the Speaking Valve (Figure 5c).

**Note:** If the Connection Strap is pulled with a certain force, the string will loosen from the clip. This is a security feature to prevent injury if the string accidentally is caught up somewhere. If this occurs, just insert the string again in the clip (Figure 5a and 5b).

**CAUTION:** Do not attempt to glue or otherwise permanently attach the string to the clip or Speaking Valve.

## Assemble speaking device and HME

**CAUTION:** Always assemble the HME with the speaking device (Speaking Valve, HME DigiTop, or HME DigiTop O<sub>2</sub>) before use (Figure 6a, 6b, 6c, 6d, 7a, and 7b). If the device is not properly assembled, the ProTrach HME 22 can get stuck in the tracheostomy tube and be difficult to remove in case of blockage.

1. Place the HME on a flat surface:  
For ProTrach HME 15 with the narrow part down and the open foam surface up (Figure 6a or 6c).  
For HME 22 with the plastic grid down and the open foam surface up (Figure 6b or 6d).
2. Put the speaking device on top of the HME and push it down until it clicks in place.
3. Check that the speaking device is firmly attached to the HME.

## 2.2 Insertion and usage of Speaking Valve with HME 15 / 22

If the Connection Strap is used, attach the clip to the neckband (Figure 8).

Attach the combined Speaking Valve and HME to the connector of the tracheostomy tube (Figure 9). **Caution:** HME 15 users should always attach the combined Speaking Valve and HME 15 to the connector gently, using a 1/4 twist in a clockwise motion. If the device is pushed too hard it can get stuck on the tracheostomy tube and be difficult to remove in case of blockage.

If you have not used HMEs previously, you should be aware that the device may increase breathing resistance to some extent.

You can now switch between speaking mode and HME mode based on your speaking needs, humidification needs, and breathing resistance needs (Figure 10). **Caution:** When changing modes, make sure that you only rotate the lid of the Speaking Valve, not the entire device. Rotation of the whole device might cause uncomfortable movement of the tracheostomy tube.

In case you feel increased breathing resistance after coughing, take out the device and wipe off the mucus, or replace the HME with a new one. Also clear the tracheostomy tube from mucus when necessary. If you feel the need to cough heavily, it is recommended that you try to remove the device from the tube prior to coughing. This prevents the device from falling

off and prevents the HME from getting clogged by mucus. If the membrane pops out through the holes in the lid of the Speaking Valve when coughing, rotate the lid, and the membrane will get back to the correct position.

## **Speaking mode**

The Speaking Valve is put in speaking mode by rotating the lid until it clicks into the speaking mode position (the flexible membrane covers the openings of the lid). In speaking mode the membrane only opens when you breathe in. It is closed when you breathe out, so you will breathe out through your mouth and nose, and you will be able to speak. You may feel a somewhat higher breathing resistance in speaking mode than in HME mode due to the membrane. You can now try to speak some words, beginning with low volume and pressure.

## **HME mode**

The Speaking Valve is put in HME mode by rotating the lid until it clicks into the HME mode position (i.e., the flexible membrane does NOT cover the openings of the lid). HME mode allows you to both breathe in and out through the HME, thereby conditioning the inhaled air.

Note that increased mucus production may appear during the first weeks of HME use.

With Speaking Valve Blue (REF 7755), the colored membrane makes it easier to see whether the valve is in speaking mode or in HME mode.

## **2.3 Insertion and usage of HME DigiTop and HME DigiTop O<sub>2</sub>**

Assemble and insert the speaking device as described in 2.1.

If you use the HME DigiTop O<sub>2</sub> always assemble and disassemble the oxygen tubing while the device is disconnected from the trach-tube. This is to avoid uncomfortable movement of the trach-tube.

Make sure that the oxygen tubing is securely attached; it should be assembled beyond the conical edge (Figure 15).

If you have not used HMEs previously, you should be aware that the device may increase breathing resistance to some extent.

When you occlude the two holes on the sides with your fingers, exhaled air will be redirected through your upper airways and you will be able to speak (Figure 11). Releasing the finger occlusion will allow you to inhale again through the

device. The two openings in the HME DigiTop / HME DigiTop O<sub>2</sub> allow you to breathe in and out through the device, thereby conditioning the inhaled air continuously.

Note that increased mucus production may appear during the first weeks of HME use.

## 2.4 Removal of the HME from the speaking devices

If you use the HME DigiTop O<sub>2</sub> remove the oxygen tubing as described above **before** you start to remove the HME. The HME is intended for single use, and must be replaced at least every 24 hours, or more often if needed. It should be removed before cleaning the speaking devices and should also be replaced if it becomes clogged by mucus.

**CAUTION:** Do not wash the HME since the important HME function then is lost.

### HME 15

1. Remove the combined speaking device and HME 15 (without separating them) from the tracheostomy tube with a gentle pulling twist in a clockwise motion.
2. Hold the speaking device with one hand, and pull the Removal Aid over the HME with the other hand (Figure 12a).
3. Pull the speaking device apart from the HME with a steady grip (Figure 12b or 12c).
4. Discard the HME 15 (Figure 12d).

### HME 22

1. Remove the combined speaking device and HME 22 (without separating them) from the tracheostomy tube.
2. Hold the speaking device with one hand, and pull the Removal Aid over the HME with the other hand (Figure 13a)
3. Press the handles together (Figure 13b) and the HME 22 will crack and loosen from the speaking device (Figure 13c or 13d).
4. Discard the HME 22 (Figure 13e).

## 2.5 Pediatric Use

If selected properly, ProTrach DualCare can be used by children of all ages, provided that it fits the local anatomy. This can change as the child grows. In children with a short neck or a high tracheostomy site, the device may contact the chin which

may be uncomfortable and can cause the Speaking Valve to twist modes inadvertently. A tracheostomy tube extender may prevent this from occurring. The ability to handle the device independently depends on the dexterity and developmental maturity of the child. Until the child is able to handle the device independently, supervision by a parent or caregiver while using the device is obligatory (see Contraindications). The indications and contra-indications for the use of the HME DigiTop and the HME DigiTop O<sub>2</sub> are the same for adults and children. However, manual occlusion for speaking may not be possible if the fingers of the child are too small. In such a case, the supervising parent or caregiver can occlude the device.

The indications and contra-indications for the use of the Speaking Valve are the same for adults and children. The ability to connect and disconnect the device from the tube, assemble and disassemble the device, and twist the device between modes depends on the dexterity and developmental maturity of the child. These tasks need to be performed by the supervising parent or caregiver until the child is able to handle the device independently.

## 2.6 Cleaning and disinfection

### 2.6.1 Cleaning of the speaking devices

In order to keep the speaking devices (Speaking Valve and HME DigiTop) clean and in a good working condition, they should be cleaned daily. If used, the Connection Strap can be cleaned together with the Speaking Valve.

**CAUTION:** Do NOT boil the device and do NOT use hydrogen peroxide for disinfection since this will damage the device.

1. Remove the HME from the speaking device as described above.
2. Flush the oxygen port connector of HME DigiTop O<sub>2</sub> (if used), with oxygen or air to remove any debris.
3. Rinse the speaking device carefully on both sides under lukewarm (20-40 °C / 68-104 °F) drinking water. Rotate the Speaking Valve lid back and forth during the rinse. Rinse the speaking device for 2 minutes (Figure 14a).
4. Mix 2 drops of dish soap in 250 ml warm (35-45 °C / 95-113 °F) drinking water. Do not use boiling water since this may damage the device.
5. Place the speaking device in the mixture for 15 minutes (Figure 14b). Make sure that the Speaking Valve is in

speaking mode (with the membrane covering the openings of the lid) before submerging it.

6. Rinse the speaking device on both sides in warm (35-45 °C / 95-113 °F) drinking water for 15 seconds. In order to reach all parts of the Speaking Valve, be sure that it is opened and closed a couple of times during rinsing (Figure 14c).
7. Inspect the device with respect to mucus and debris residues. Repeat the steps above if necessary.
8. Flush the oxygen port connector of HME DigiTop O<sub>2</sub> with oxygen or air to remove any remaining water.
9. After cleaning, put the speaking device on a clean gauze and leave it to air-dry for at least 2 hours (Figure 14d). Make sure that the Speaking Valve is in speaking mode.
10. When the speaking device is completely dry, it can either be used again or placed in a clean, closed container for storage.

### 2.6.2 Disinfection of the speaking devices

It is recommended that the device is disinfected at least once a week or before use if it looks dirty or there has been a risk of contamination. A risk of contamination could exist if the device has been dropped on the floor or if it has been in contact with a pet, someone with a respiratory infection, or any other gross contamination. If disinfection is needed always clean the device first as described above and then disinfect by using the following procedure:

1. First follow the cleaning procedure as described in 2.6.1, without step 9 and 10. After cleaning, put the Speaking Valve in speaking mode and submerge it for 10 minutes either in ethanol 70% or isopropyl alcohol 70% (Figure 15a).

**CAUTION:** Do NOT use hydrogen peroxide.

2. After cleaning and disinfection, put the speaking device on a clean gauze and leave it to air-dry for at least 2 hours (Figure 15b). Make sure that the Speaking Valve is in speaking mode.

3. When the speaking device is completely dry, it can either be used again or placed in a clean, closed container for storage.

**CAUTION:** Do not use the speaking device until it is completely dry. Inhalation of disinfectant fumes can cause severe coughing and irritation of the airways. The HME DigiTop O<sub>2</sub> oxygen port connector shall be flushed with oxygen or air before use to remove any disinfection residuals.

## 2.7 Storage instructions

When not in use, clean and disinfect the Speaking Valve as described above, and then store it in a clean and dry container at room temperature. Protect from direct sunlight.

## 2.8 Device lifetime

The HME is for single use and must be replaced at least every 24 hours, or more often if needed.

The Speaking Valve, HME DigiTop, HME DigiTop O<sub>2</sub>, and Connection Strap may be used for a maximum of 2 months as long as they are intact and are functioning as intended.

## 2.9 Disposal

Always follow medical practice and national requirements regarding biohazard when disposing of a used medical device.

## 2.10 Accessories

**ProTrach Connection Strap:** Can be used to attach the Speaking Valve to the neckband.

**Removal Aid:** Can be used to help to remove the HME from the speaking device.

## 3. Additional information

### 3.1 Compatibility with MRI Examination

**MR-Safe:** This device does not contain any metallic elements and has no potential for interaction with the MRI field.

### 3.2 Ordering information

See end of this Instructions for Use.

### 3.3 User assistance information

For additional help or information, please see back cover of this Instructions for Use for contact information.

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## Ordering information

Name	REF	Rx
ProTrach DualCare Set 15: 30 pcs HME 15 Regular+ 1 pc Speaking Valve + 1 pc HME DigiTop + 1 pc Removal Aid + 1 pc Connection Strap	7741	Rx
ProTrach DualCare Set 22: 30 pcs HME 22 Regular+ 1 pc Speaking Valve + 1 pc HME DigiTop + 1 pc Removal Aid + 1 pc Connection Strap	7740*	Rx
ProTrach DualCare Speaking Valve: 1 pc Speaking Valve + 1 pc HME DigiTop + 1 pc Connection Strap	7744	Rx
ProTrach DualCare Speaking Valve Blue: 1 pc Speaking Valve Blue + 1 pc HME DigiTop Blue + 1 pc Connection Strap	7755	Rx
ProTrach HME 15 Regular, 30 pcs	7742	Rx
ProTrach HME 15 XtraMoist, 30 pcs	7766	Rx
ProTrach HME 22 Regular, 30 pcs	7747*	Rx
HME DigiTop O2, 1 pc	7756	Rx
Removal Aid, 2 pcs	7745	
ProTrach Connection Strap, 2 pcs	7746	

\* not available in all countries

**Caution:** United States Federal law restricts the sale, distribution and use by, or by order of a physician or a licensed practitioner of devices identified herein with the symbol “Rx”. The prescription-free availability of these products outside the United States may vary from country to country.

Unregistered copy

Atos Medical AB  
Kraftgatan 8, P.O. Box 183, SE-242 22 Hörby, Sweden  
Tel: +46 (0)415 198 00 • Fax: +46 (0)415 198 98  
[www.atosmedical.com](http://www.atosmedical.com) • [info@atosmedical.com](mailto:info@atosmedical.com)

© Atos Medical AB, Sweden.