FreeHands FlexiVoice™

Instructions for Use

Rx
ONLY

Atos

Atos Medical Your voice


Figure 1


Figure 2


Figure 3


Figure 4


Figure 5


Figure 6


Figure 7


Figure 8


Figure 9


Figure 10


Figure 11


Figure 12


Figure 13


Figure 14


Figure 15


Figure 16


Figure 17


Figure 18


Figure 19


Figure 20


Figure 21


Figure 22


Figure 23

Prescription information

CAUTION: United States Federal law restricts this device to sale, distribution and use by or on order of a physician or a licensed practitioner. The availability of this product without prescription outside the United States may vary from country to country.

Disclaimer

Atos Medical offers no warranty – neither expressed nor implied – to the purchaser hereunder as to the lifetime of the product delivered, which may vary with individual use and biological conditions. Furthermore, Atos Medical offers no warranty of merchantability or fitness of the product for any particular purpose.

Patents and trademarks

Provox® is a registered trademark owned by Atos Medical AB, Sweden.

FreeHands FlexiVoice™ is a trademark of Atos Medical AB.

For information about protective rights (e.g. patents), please refer to our web page www.atosmedical.com/patents.

Unregistered copy

Contents

| | |
|----------------------------|-----|
| ENGLISH | 7 |
| DEUTSCH | 13 |
| NEDERLANDS | 20 |
| FRANÇAIS | 27 |
| ITALIANO | 34 |
| ESPAÑOL | 41 |
| PORTUGUÊS | 48 |
| SVENSKA | 55 |
| DANSK | 61 |
| NORSK | 67 |
| SUOMI | 73 |
| MAGYAR | 79 |
| EESTI | 85 |
| HRVATSKI | 91 |
| РУССКИЙ | 97 |
| 日本語 | 104 |
| 한국어 | 110 |
| 简体中文 | 116 |
| 繁體中文 | 121 |
| العربية | 131 |
| ORDERING INFORMATION | 132 |

The Instructions for Use, which accompanies this product, may be revised from time to time and must therefore be reviewed prior to each procedure in which the product is used.

1. Descriptive information

1.1 Intended use

The Provox FreeHands FlexiVoice combines pulmonary rehabilitation using a Heat and Moisture Exchanger with voice rehabilitation using an Automatic Speaking Valve or Manual Occlusion, in laryngectomized patients using a voice prosthesis.

1.2 CONTRAINDICATIONS

The Provox FreeHands FlexiVoice is not intended to be used by patients unable to remove or operate the device, unless the patient is under constant supervision of a clinician or a trained caregiver. For example, patients who are unable to move their arms, patients with decreased levels of consciousness, or patients with diseases that put them at risk for unpredictable periodic loss of consciousness.

1.3 Description of the devices

The Provox FreeHands FlexiVoice consists of two parts assembled together, a speaking valve for single patient use and a disposable HME cassette (Fig. 1).

The speaking valve is made of plastic and the membrane is made of silicone. The HME cassette is also made of plastic and a salt treated polyurethane foam.

The speaking valve has two modes; Automatic Speaking Mode and Locked Mode (Fig. 2). Rotating the top of the speaking valve moves the device into the automatic speaking or the locked position. Speaking can be done both by using the automatic speaking valve and by manual occlusion of the opening in the front. Manual occlusion is possible in both modes.

The speaking valve is in automatic speaking mode when the flexible membrane is not locked by the hook (Fig. 3). In automatic speaking mode, the membrane can be closed automatically by the exhaled air. The air is then redirected through the voice prosthesis. The membrane opens again automatically when breathing in.

In locked mode, the flexible membrane is locked by a hook that goes into a loop in the membrane (Fig. 4). This prevents the membrane from closing when breathing heavily, for example, during physical activity.

There are four different speaking valves, each with a different strength membrane: light, medium, strong and xtrastrong.

On the membrane part that protrudes on the top of the device, there are 1, 2, 3 or 4 dots visible. The number of dots represents the strength of the membrane (Fig. 5).

1 dot = Light membrane (easiest to close)

2 dots = Medium membrane

3 dots = Strong membrane

4 dots = XtraStrong membrane

The Arch is an accessory that can be used with the Provox FreeHands FlexiVoice to keep fabric away from the opening of the device. It is simply pressed onto the front of the Provox FreeHands FlexiVoice (Fig. 6). Please note that in some positions the Arch may restrict the possibility of manual occlusion.

The Removal Aid (Fig. 7) can be used to remove the HME cassette from the speaking valve.

Content of the sets

Provox FreeHands FlexiVoice Set & Provox FreeHands FlexiVoice Set Plus

- 3 Speaking Valves (Light, Medium, and Strong)
- 1 Removal Aid
- 2 Arches
- 30 Provox FreeHands HME cassettes
(Provox FreeHands FlexiVoice Set Plus only)

Provox FreeHands FlexiVoice Light, Medium, Strong and XtraStrong

- 1 Speaking Valve (of the chosen strength)
- 1 Arch
- 1 Removal Aid

1.4 WARNINGS

- Do not use the HME cassette without assembling it first with Provox FreeHands FlexiVoice Speaking Valve or Provox HME Cap. The HME can get stuck in the attachment holder when used by itself (Fig. 8).
- Do not use Provox FreeHands FlexiVoice with any other attachment devices than original Provox HME attachment devices. The product can get stuck in the attachment device. The attachment device then has to be removed or may be damaged, possibly leading to restricted breathing.
- Do not reuse any part of the Provox FreeHands FlexiVoice between patients. This may cause cross contamination. The device is intended for single patient use only.

1.5 PRECAUTIONS

- Do not modify the product. The speaking function and/or HME function will be lost.
- Always clean the speaking valve according to the instructions below. If the product is not properly cleaned, it may cause infection and product malfunction.
- Do not administer medicated nebulizer treatment over the device since the medication can become deposited in the device.
- Do not use humidifiers or heated humidified oxygen via a mask over the tracheostoma while using the device. The HME will become too wet. If oxygen therapy is required, use only non-heated humidified oxygen.

2. Instructions for use

2.1 Checkout procedures

Before each use, check that the speaking valve and HME cassette are undamaged and that the speaking valve functions as intended. It should easily turn into the locked mode and automatic speaking mode (Fig. 2). If the device does not function as intended or seems damaged, do not use the device and get a replacement.

2.2 Operating instruction

Do not use the speaking valve during sleep. The opening of the speaking valve may be covered or the speaking valve may unintentionally be turned to speaking mode resulting in restricted breathing. During sleep, use a Provox XtraHME cassette or a Provox HME Cap assembled with a HME cassette (Fig. 9).

Membrane strength selection

Try the Medium membrane (2 dots) first. If you experience unintended closing of this membrane during normal breathing or light activity, try the Strong or the XtraStrong membrane (3 or 4 dots). If you feel it takes too much effort to close the Medium membrane for speaking, try the Light membrane (1 dot). Make sure the membrane does not close unintentionally during regular breathing.

You can also use a stronger membrane when you are more active and a lighter one when you are less active.

Assembly

Place the HME cassette on a flat surface with the plastic grid down and the open foam surface up (Fig. 10). Place the speaking valve on top of the HME cassette and push it straight down until it clicks in place. Check that the speaking valve is firmly attached to the HME cassette.

Caution: Use of non-Provox accessories may cause product malfunction and harm to the user.

Insertion and usage

Insert the Provox FreeHands FlexiVoice into the connector of the Provox HME attachment device (Fig. 11). Breathe normally. Switch the device between the modes by rotating the top of the speaking valve until it stops in each position (Fig. 12).

Caution: Rotate only the top of the speaking valve when changing modes, not the entire device. Rotating the entire device can cause uncomfortable movement of the attachment device.

Automatic Speaking Mode

The Provox FreeHands FlexiVoice is in Automatic Speaking Mode when the flexible membrane can move freely and is not locked by the hook (Fig. 3).

In this mode, the increase in exhaled airflow for speaking automatically closes the membrane.

Manual Occlusion

The Provox FreeHands FlexiVoice can also be manually occluded for speaking by occluding the front opening with a finger (Fig. 13). To breathe in, release the finger occlusion so that the front opening is open again. Manual occlusion can be used in both Automatic Speaking Mode and in Locked Mode.

Coughing

In automatic speaking mode, the membrane can pop out when coughing. If this happens, breathe normally and push the membrane back through the opening with your fingertip (Fig. 14). If necessary, remove the Provox FreeHands FlexiVoice and re-position the membrane.

Removal of the HME cassette from the Speaking Valve

The HME cassette is intended for single use. It should be removed at the end of each day before cleaning of the speaking valve. It should also be replaced if it is clogged by mucus and the breathing resistance has increased.

Caution: Do not attempt to wash the HME cassette since this will wash away the impregnation essential for its function.

1. Hold the speaking valve (or cap, if used) with one hand, place the Removal Aid over the HME cassette with the other hand (Fig. 15).
2. Press the handles together (Fig. 16). The HME cassette will crack and loosen (Fig. 17).
3. Discard the HME cassette (Fig. 18).

Cleaning

In order to keep the speaking valve in a good working condition and to avoid build-up of microorganisms, it should be cleaned and disinfected daily.

Caution: Do not clean or disinfect the speaking valve in any other way than described below. Do NOT boil the device and do NOT use hydrogen peroxide for disinfection since these actions will damage the device.

1. Remove the HME cassette from the speaking valve as described above.
2. Rinse the speaking valve carefully on both sides under lukewarm (20-40 °C / 68-104 °F) drinking water while rotating the lid back and forth for 2 minutes (Fig. 19).
3. Mix 2 drops of dish soap in 250ml warm (35-45 °C / 95-113 °F) drinking water. Do not use boiling water since this may damage the speaking valve.
4. Put the speaking valve in speaking mode and place the speaking valve in the mixture for 15 minutes (Fig. 20).
5. Rinse the speaking valve on both sides in warm (35-45 °C / 95-113 °F) drinking water for 15 seconds (Fig. 21). In order to reach all parts of the speaking valve, be sure that it is switched between speaking and locked mode a couple of times during rinsing.
6. Inspect the device regarding mucus and debris residues. Repeat the steps above if necessary.
7. After cleaning, place the speaking valve in speaking mode on clean gauze and leave it to air-dry for at least 2 hours (Fig. 23). When the speaking valve is completely dry, it can either be used again or be placed in a clean, closed container for storage.

Caution: Do not use the speaking valve until it is completely dry. Inhalation of disinfectant fumes can cause severe coughing and irritation of the airways.

Disinfection of the speaking valve

Disinfection might be needed if the speaking device looks dirty or there has been a risk of contamination. A risk of contamination could exist if the device has been dropped on the floor or if it has been in contact with a pet, someone with a respiratory infection, or any other gross contamination. Then the device should be both cleaned and disinfected before use. If disinfection is needed first clean the device as described above and then disinfect by using the following procedure:

1. First follow the cleaning procedure as described in section "Cleaning" above, without step 7. After cleaning, put the speaking valve in speaking mode and submerge it for 10 minutes either in ethanol 70% or isopropyl alcohol 70% (Fig. 22).

Caution: Do NOT use hydrogen peroxide.

2. After cleaning and disinfection, put the speaking valve on a clean gauze and leave it to air-dry for at least 2 hours (Fig. 23). Make sure that the speaking valve is in speaking mode.
3. When the speaking valve is completely dry, it can either be used again or placed in a clean, closed container for storage.

Caution: Do not use the speaking valve until it is completely dry. Inhalation of disinfectant fumes can cause severe coughing and irritation of the airways.

2.3 Storage instructions

When not in use, clean and disinfect the speaking valve as described above and then store it in a clean and dry container at room temperature. Protect from direct sunlight.

2.4 Device lifetime

The Provox FreeHands FlexiVoice Speaking Valve can be used for a maximum of 6 months as long as it is intact and functions as intended. Usage beyond this time may lead to improper valve function, e.g., delayed closing and/or opening, air leaks and increased closing sound, and/or other changes in the function due to device ageing.

The HME cassette is for single use.

2.5 Disposal

Always follow medical practice and national requirements regarding biohazards when disposing of a used medical device.

2.6 Accessories

Caution: Use only original components. Other accessories may cause product malfunction.

Provox HME Cap is a reusable dome-shaped ring that allows you to use HME cassette during sleep, or without the use of the Provox FreeHands FlexiVoice Speaking Valve. Provox HME Cap is made of titanium.

Provox HME attachment devices are products that are used to hold products of the Provox HME system in front of the stoma, e.g., Provox Adhesives, Provox LaryTubes and Provox LaryButtons.

Removal Aid can be used to help to remove the HME from the speaking valve.

3. Troubleshooting information

In case of problems during use, remove and inspect the device. Below are some common problems and the recommended actions.

“The device leaks” – Inspect the speaking valve for leakage or damage. Replace the device if it is damaged. Clean the device as described in section “Cleaning and disinfection”. Also check the attachment device and HME cassette and replace if needed. If there is still leakage after the cleaning process, replace with a new speaking valve.

“The device closes if I increase breathing just slightly” – Use a device with a stronger membrane. Alternatively, use a stronger membrane when physical activity is expected.

“The membrane closes when I breathe heavily” – Unintended closing can occur while breathing heavily. To prevent this, use the locked mode position for the device in foreseeable situations. Speaking is still possible using manual occlusion. Unlock the device again to enable handsfree speech (Fig. 2).

“The membrane doesn’t close” – First check that the device is set in speaking mode. Inspect the speaking valve and membrane position. The membrane could be stuck behind the hook. This may occur after being re-positioned wrongly when coughed out. Gently position the membrane correctly between the hook and hole (Fig. 3). If you find any damage to the speaking valve, replace it.

“The membrane pops out when I speak” – The membrane is too light for you. Try using a stronger membrane.

“The membrane pops out when I cough” – This is normal. Remove the device with the attached HME cassette from the adhesive. Put the membrane back in position with your finger (Fig. 14). If this occurs frequently, you may want to try a stronger membrane.

“The speaking valve components came apart” – Do not try to re-assemble. Replace the device.

“It is hard to breathe” – The HME cassette may be clogged by mucus. Change to a new HME cassette and check again. If the problem persists, clean the device according to section “Cleaning and disinfection”. Inspect the device and make sure it is undamaged. Replace if damaged.

“The clothes I use to cover my neck interfere with my breathing” – Assemble the provided Arch to the speaking valve (Fig. 6).

“I hear a popping noise when I start to speak” – Try to initiate your speech in a gentler, less abrupt manner. Try a stronger membrane that closes less abruptly.

“I hear air escaping when I start to speak” – Try to initiate your speech in a more sudden way so that the membrane closes faster. Try a lighter membrane that closes faster.

4. Additional information

4.1 Compatibility with MRI Examination

MR-Safe: This device does not contain any metallic elements and has no potential for interaction with the MRI field.

4.2 Travel or international use

Provox FreeHands FlexiVoice and the HME are metal free and x-ray translucent. When traveling to other countries, ensure that you have products with you. Also ensure you clean the device in drinking water. Visit www.atosmedical.com and/or contact your local Atos Medical office for contact information for replacement products.

4.3 User assistance information

For additional help or information, please see back cover of the Instructions for Use for contact information.

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© Atos Medical AB, 2018.