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📦 Packaging Unit

Purpose:

The maintenance tubes are intended for use in patients with dysphagia. They are comprised of an outer tube and two inner cannulae – some with universal attach-ment, cuff for blocking and/or access for suction cleaning. The maintenance cannulas are suitable for use during radiotherapy and MRI ex-aminations.

Description:

Prior to insertion of the tube the cuff must be inflated using a cuff-pressure measur-ing device. If the cuff is securely tight it is re-moved; the maintenance tube can then be inserted by a physician or nurse or other trained person. To ensure easy insertion of the maintenance tube into the tracheos-tomy, it should be treated with a specially-designed lubricant. When blocking the cuff, make sure not to use pressure greater than 25 mmHg (equivalent to 34 cmH₂O) over a longer period. That is why the cuff pressure should be checked at regular intervals using a cuff pressure measuring de-vice. If the tube becomes clogged due to pronounced secretion, the inner cannula must be removed and cleaned. Then use a tracheal suction catheter to suction-clean it through the outer tube. The inner cannula is then re-inserted.

When removing the maintenance tube, any fluid which has gathered above the cuff must be suctioned off before the cuff is released.

The upper section of the inner cannulae is coloured. Based on the colour code, the inner cannula can be matched with the fitting outer tube.

The inner cannula and outer tube must be matched as follows:

Size	Inner dia. outer tube	Outer dia. inner cannula	Colour code inner cannula	Length
7	7 mm	6.8 mm	green	80 mm
8	8 mm	7.8 mm	white	80 mm
9	9 mm	8.8 mm	yellow	85 mm
10	10 mm	9.8 mm	blue	95 mm
11	11 mm	10.8 mm	pink	97 mm

Maintenance tubes are single-patient products. They are made of medical grade plastic material. The maximum period of use is 29 days. The maintenance tube should be removed after 6 days of wearing time at the most for a tracheal inspec-tion, for example for prevention and early recognition of stenoses. Before re-inser-tion the maintenance tube should be cleaned as described in the section “Cleaning”.

Cleaning:

The maintenance tube is cleaned without chemical additives in a saline solution at a maximum of 65°C. Following this we recommend disinfection of the inner can-

nulae and outer tubes using a suitable disinfectant following the instructions for use. The inner cannula is then treated for maintenance with a lubricant (e.g. **PRIMA-STOM II Oil** using a soft cloth and following the instructions for use. Follow-ing cleaning and disinfection, the cuff must be tested to make sure it does not leak.

Safety information:

The product must not be used if the packaging is damaged. Non-compliance with the cleaning instruction may lead to candida contamination because the inner can-nula was no longer seated tightly in the outer tube. Excessively high blocking pres-sure in the cuff may cause pressure necroses in the trachea.

Release the maintenance tube completely before introducing it into the trachea. Prior to blocking, the maintenance tube should not be retracted to prevent obstruc-tion of the inner diameter by the block.

Caution: Leaving the protective cap on the pilot balloon during use may result in a loss of cuff pressure.

The protective cap can be discarded once the tube has been unpacked. It is only used for initial sterilisation by the manufacturer.

The product is not suitable for use with laser devices. Serious incidents involving the single-use inner cannulas must be reported to the manufacturer and the competent authority.

Sterility information:

The products have been sterilised using ethylene oxide. Sterility is guaranteed only if the packaging is unopened, undamaged or torn. The products may not be re-sterilised.

Product disposal:

Disposal of maintenance tubes and their accessories must be carried out in accord-ance with current national regulations applying to waste products.

Caution: Warranty, guarantee and liability claims against the manufacturer will be nullified if any changes are made to the product by users or if the product is used outside of the scope of applications described in these instructions.

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