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 Packaging unit

Purpose:

HME (Heat and Moisture Exchangers) moisturise, warm and clean the respiratory air of laryngectomy and tracheostomy patients. In this way, the formation of stubborn mucus in the trachea and the risk of infection are minimised. Additional oxygen can be administered to the patient using the oxygen attachment.

indications:

Patients with spontaneous breathing via a tracheal tube, in the hospital or a domestic environment.

Instructions:

Remove the HME from its packaging and connect the round opening to the attachment ring of the cannula or of the plaster or to the universal connector (illustrations 1 and 2). When breathing out, the HME assumes the moisture and warmth of the exhaled air and re-introduces these into the inhaled air. This moisturising and warming effect activates after a few breaths.

Contra-indications:

The HME is not to be used outside the tidal volume range (50 – 1000 ml.) as an increased dead-space in strongly reduced tidal volumes results in CO₂ retention or, in strongly increased tidal volume, in insufficient moisturisation. Usage in patients with significantly reduced fluid intake (exsiccosis) is not permitted, nor in cases of excessive secretions in the lung or airways or with a high level of air leakage (expiratory air does not flow through the HME.)

Cleaning:

The HME is a disposable product and, for hygiene reasons, should be replaced daily. If respiratory resistance is increased due to clogging of the HME with secretions, replacement must be carried out earlier. For hygiene reasons, cleaning of the HME is not recommended. Re-sterilisation of the sterile product version is not permitted, as this can affect its functioning.

Safety Information:

Prior to device operation, check the connection between the HME and the corresponding attachment to make sure it is firmly seated and leak-tight. If the HME is blocked or shows other signs of malfunction it must be removed immediately and disposed of. It must not be re-used or cleaned, and the use of disinfectants can render an HME unusable. A risk of infection can arise if the HME is used for too long or if unauthorised cleaning procedures are performed. The patient should be supervised at regular intervals for additional respiration needs (indication: multiple thick secretions). The HME must not be used together with warm-water moisturisers or nebulisers.

Storage:

HME can be stored under normal ambient conditions (5 - 35 °C. and 30 - 70 % relative humidity).

Product disposal:

Disposal of products must be carried out in compliance with current national regulations applying to waste products.

Caution: Warranty, guarantee and liability claims against the manufacture are nullified if changes are made to the product by users or if the product is used outside of the scope of applications described in these instructions.

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