PRIMA-DYS®

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UK Instructions for use

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 Packaging unit

Purpose

The **PRIMA-DYS**[®] dysphagia tube is specially designed for patients with swallowing difficulties of all aetiologies and severity in protecting the lower airways. Among these may be included ventilated patients in weaning (the process of withdrawal from ventilation).

Optimal tracheal hygiene and the resultant improvement in pneumonia prophylaxis is obtained using the integral sub-glottal suction system.

Description

Prior to insertion of the tube the cuff must be inflated using a cuff-pressure measuring device to check for tightness. If securely tight the cuff may be deflated and the

PRIMA-DYS[®] inserted by medical staff. To ensure easy in-tertion into the tracheostomy the **PRIMA-DYS**[®] should be treated with a specially-de-gined lubricant. When inflating the cuff, care must be taken to ensure the pressure never exceeds 25 mm. Hg for long periods. The cuff-pressure should therefore be checked regularly using a cuff-pressure measuring device. When removing the **PRIMA-DYS**[®] any fluid which has gathered above the cuff must first be suctioned before the cuff is released. By means of a suction opening directly above the cuff, any aspirated, accumulated secretions can be regularly removed through the outlet suction-tube. Here the cuff remains filled with air (sub-glottal suctioning in the case of a blocked tracheal cannula), thus optimising protection of the lower airways. In the case of a blocked tracheal cannula suctioning via the sub-glottal suction system can be achieved using a syringe. This process facilitates better documentation of the aspirated secretion volume in a certain time-frame and the identification of the suctioned material by appearance and type (purulent, clear, bloody, mixed with food or tube-fed food, etc.).

Use of the suctioning protocol can be helpful here, finding the best time for a safe unblocking with the highest possible safety of the lower airways. For improved tracheal hygiene (e.g. softening of hardened surfaces, dilution of thick secretions - or similar) and under the doctor's directions/instructions sub-glottal flushing with a sterile saline solution (0.9% NaCl) can be carried out.

Before such a sub-glottal flushing the cuff-pressure is again checked using a cuff-pressure measuring device to obtain the best pressure in protecting the lower airways against aspiration. Any accumulated secretions above the cuff must be suctioned away. Using a syringe the fluid is directed over the suction tube and immediately suctioned away again with a syringe.

In order to keep the risk of pneumonia as low as possible the **PRIMA-DYS**[®] functions without an inner cannula (a tightness feature in cases of chronic aspiration of saliva or secretion); attempts at deflating and speaking are therefore possible.

By means of the sub-glottal suction system the **PRIMA-DYS**[®] facilitates deflation of the cuff, ensuring the best possible protection of the lower airways.

Prior to release of air from the cuff, accumulated secretions over the cuff are removed using a syringe and preferably documented – as described above.

After no more secretions can be removed, air is released from the cuff, during which the trachea is suctioned in accordance with established procedure. In this way, any migration into the lower airways of incompletely removed secretions above the cuff – or sudden aspirated material as a result of reflex or vomiting – is prevented.

PRIMA-DYS[®] is a one patient only produkt. The cannula adapts to the anatomy of the trachea by virtue of the plastic material used. Maximum time in situ is 29 days. After a maximum time in situ of 6 days the **PRIMA-DYS**[®] should be removed for inspection so that e.g. stenosis may be avoided or recognised early. Before re-insertion the cannula should be cleaned as described in Section 1.3.

Caution: The product is not designed for use with lasers.

Cleaning

The **PRIMA-DYS**[®] is cleaned without chemical additives in a saline solution (0.9% NaCl), at a maximum temperature of 65 °C. We recommend that this be followed by disinfection of the tube using a dedicated solution in accordance with the usage instructions. After successful cleaning and disinfection the cuff must be tested for tightness. Failure to adhere to the cleaning instructions invalidates the product guarantee.

Safety information

If the packaging is damaged the product must not be used. If the cleaning instructions are not followed, a build-up of candida deposits may occur. If the cuff is positioned too tightly, pressure necrosis in the trachea can result. Completely unblock the **PRIMA-DYS**[®] before introducing it into the trachea.

Prior to blocking, the **PRIMA-DYS**[®] should not be retracted to prevent obstruction of the inner diameter by the block.

If the cuff is positioned too tightly, pressure necrosis in the trachea can result.

Serious incidents involving the **PRIMA-DYS**[®] dysphagia tube must be reported to the manufacturer as well as to the competent authority.

Sterility information

The products have been sterilised using ethylene oxide. Sterility is guaranteed only if the packaging is unopened, undamaged or torn. The products may not be re-sterilised. Proceed with special care when selecting and placing the **PRIMA-DYS**[®] to avoid complications by forming granulation tissue.

Product Disposal:

Dispose of the product must be carried out in compliance with current national regulations applying to waste products.

Caution: Warranty, guarantee and liability claims against the manufacturer will be nullified if any changes are made to the product by users, if the product is used outside of the scope of applications described in these instructions for use or if the cleaning instructions are not complied with.