


Tracoe Percutaneous Dilatation – REF 520

Tracoe Experc Set - For percutaneous dilation based on Ciaglia technique.

Product Description

Percutaneous Dilatation Tracheostomy (PDT) has become routine in intensive care. The Ciaglia technique is now the method most widely used. The Tracoe Experc system is based on this proven technique. The products are tailored to the use of fenestrated and non-fenestrated tracheostomy tubes, including tubes with subglottic suction. They are supplied already mounted on their minimally traumatic insertion systems (P-tubes).



General Information

| | |
|--------------------|--|
| Product line | Percutaneous Dilatation |
| Legal manufacturer | Tracoe Medical GmbH Reichelsheimer Str. 1/3 55268 Nieder-Olm |
| Country of origin | Germany |
| Handling | |
| Shelf life | 5 years, as from date of production |

Material Specifications

| | |
|---------------------------|--------------------------------|
| Material short dilatator | Polypropylene |
| Material Guiding Catheter | Fluorinated ethylene propylene |
| Material experc Dilator | Polyurethane/ PVP |
| Material guide wire | Nitinol, TFE coated |
| Biocompatibility | EN ISO 10993-1 |



Storage




Packaging

| | |
|-----------------------------|------------------------------------|
| Packaging unit | 1 Unit |
| Weight | 495 g |
| Dimensions (W x H x L) | 164 x 430 x 38 mm |
| Material outer packaging: | Cardboard box |
| Material sterile packaging: | Hard blister sealed in sterile bag |

Certification

| | |
|---------------------------------|------------------------------------|
| Product classification acc. MDR | IIa |
| CE Marking | |
| GMDN | 14099 |
| Quality system certified to | EN ISO 13485 |
| FDA registration | Clearance number K043160 / D020866 |
| Global Trade Item Number | 04035324017263 |



Scope of Supply

- 1 Scalpel
- 2 Syringe
- 3 Puncture needle (14 G) with teflon catheter
- 4 Seldinger guide wire made of kink resistant Nitinol with inserter
- 5 Short dilatator (14 Fr)
- 6 Guiding catheter with safety stop
- 7 TRACOE expirc dilatator with hydrophilic coating
- 8 4 compresses



Tracoe Medical GmbH
Reichelsheimer Str. 1 / 3
55268 Nieder-Olm
Germany



Breathing-Speaking-Living

Atos Medical AB
Hyllie Boulevard 17
215 32 Malmö
Sweden
www.atosmedical.com